

Dose Calibrators

Increased Functionality for Your Dose Calibrator...

CRC[®]-15R Enhanced Software

Add more functionality to your CRC-15R Dose Calibrator with the CRC-15R ENHANCED SOFTWARE, with all of the same great features of the CRC-15R plus a few new ones.

- ★ Inventory control
- ★ Dose table calculations
- ★ Newly added linearity protocol
- ★ Radiochromatography QC package
- ★ Built in geometry testing

The CRC-15R ENHANCED SOFTWARE will still perform the Daily QC protocol, zero adjust, background subtraction, system test, accuracy (with a deviation percentage) and constancy. In addition, the software will calculate dose tables, perform geometry tests (for between 2 and 10 measurements), do radiochromatography, calculate linearity (using standard, Lineator or Calicheck methods) and maintain a radiopharmaceutical inventory of kits or unit doses. All of the acquired data are printable.

The CRC-15R Enhanced Software is an available option in all new CRC-15R Dose Calibrators or it can be added to your existing CRC-15R.

15R ENHANCED FEATURES

- ★ Capintec quality and reliability
- ★ Calculation of a printable dose table
- ★ Inventory control - maintains inventory of unit doses and kits
- ★ Linearity test - calculates a report for standard, lineator and calicheck methods
- ★ Geometry test showing variations for 2 to 10 different measurements
- ★ Radiochromatography of kit QC strips - one strip, two strip, HMPAO and MAG3

STANDARD LINEARITY		
#	MEASURED	%VAR
1	50.4mCi	2.10
2	36.5mCi	1.85

Any Key to Continue

Standard Linearity Test Results Screen

INVENTORY	
1.	ADD
2.	KIT
3.	WITHDRAW
4.	PRINT
5.	DELETE

Inventory Menu

QC TEST	
Single Strip	
Top	/ (Top + Bottom) %
Bottom	/ (Top + Bottom) %

Any Key to Continue

Single Strip Result Screen

