

# In Step With Capintec

www.capintec.com

The Newsletter of Capintec, Inc.

Winter 2006

## In This Issue

**Capintec Signs Agreement with Numa**  
...page 1

**Capintec RSNA 2006 Preview**  
...page 2

**Employee Updates**  
...page 3

**Show Schedule**  
...page 3

**Lead-Lined, Steel Frame Cabinetry**  
...page 4

## Capintec Partners with Numa for DICOM Compatibility

NumaCap™ Software to integrate the Captus® 3000 Thyroid Uptake System with RIS and PACS systems

*Capintec is proud to announce a new software application for its Captus® 3000 Thyroid Uptake System designed to meet the needs of our customers working in a digital reporting environment*

Capintec, a leading supplier of nuclear medicine and PET products and services; and Numa, Inc, the nuclear medicine connectivity company, has signed an agreement to develop an application for the digital reporting and communication of Captus® 3000 Thyroid Update Systems. These systems measure thyroid function utilizing radionuclide uptake for data.

Under the agreement, Numa will develop a digital application for integration into the thyroid uptake device that will create a complete DICOM record for each patient exam, including demographics and a secondary screen capture of the alphanumeric exam data. The new NumaCap software will interface with a RIS or PACS and query the DICOM worklist provider to populate the Captus® report with patient demographics, saving facilities significant time in data re-entry.

After study completion, the Numa application will export a snapshot of the test results. Completed patient DICOM files will be available for viewing on any DICOM nuclear medicine or PACS workstation and for transmission over the hospital network or Internet as well as archiving using a PACS.

"As medical facilities increasingly convert to a digital environment with connectivity across the enterprise, digital storage and transmission of all nuclear medicine exams is crucial," commented Lawrence Smith president of Numa. "Numa is proud to have been selected by Capintec, the

industry leader, to enable this advanced connectivity for thyroid uptake exams and looks forward to a long and productive relationship with the company."

Capintec's Captus® 3000 Thyroid Uptake System software provides fast, dedicated modules for Thyroid Uptake, Wipe Tests, Schilling Tests, Blood Volume, RBC Survival, and Bio-Assay. It comes equipped with a "Quick Start-Up" feature that walks the user through all program setup functions that not only helps the user perform patient studies quickly, but also serves as a training tool for additional users.

"As the leader in thyroid update devices, we are pleased to offer new state-of-the-art digital technology that will eliminate the manual

entry of patient data and streamline reporting and patient file management," said John Viscovic, VP Sales for Capintec. "As a leader in nuclear medicine connectivity, Numa is an excellent partner to firmly position this device in today's digital age."

### About Numa

For over fourteen years, Numa has been the leader in nuclear medicine connectivity and system integration, with solutions that deliver enhanced throughput, efficiency and productive workflow. Numa meets specific needs in nuclear medicine, nuclear cardiology, PET and PET/CT for image management, archiving, multi-vendor connectivity, DICOM modality worklist management, and secure file transfer throughout the department, enterprise and beyond. Visit their website at: [www.numa-inc.com](http://www.numa-inc.com) ■



"Quick Start-Up" Menu



### **CRC®-Ultra Dose Calibrator**

Stop by to see the CRC®-Ultra Dose Calibrator, along with a Nuclear Medicine Chamber, PET Chamber, Well and Beta Counter. Demonstrating just some of the many combinations that are possible with the CRC®-Ultra, Capintec representatives will be available to help answer your questions and help you select the right combination for your department.

While at the booth, don't forget to pick up a free CRC®-Ultra brochure, showcasing its various combinations, components and specifications.

**See These Items and More  
at RSNA Booth #8929!**

Our Lead-Lined Cabinetry features solid steel frames: an industry exclusive! Built to your specifications, Capintec can help guide you in the design process of your new or upgraded cabinetry needs. Made with the ideals of safety, ease-of-use and longevity, these durable cabinets come with a lifetime guarantee.



Stop by our booth to see a demonstration model, and pick up a copy our new brochure!

### **Lead-Lined Cabinetry**

# CAPINTEC RSNA

# 2006 PREVIEW



As the exclusive distributor of Radiology Support Devices, Inc.'s Heart, Thorax, Striatal, and Thyroid Phantoms, Capintec will be displaying the RSD Pixy™ Transparent Phantom at RSNA. Speak with a Capintec representative for more information!

### **RSD Pixy™ Transparent Phantom**

Perform Thyroid Uptake measurements quickly and easily with the Captus® 3000 - see it in person at our RSNA booth! With a user-friendly "Quick Start-Up" feature, fast thyroid uptake modules for Thyroid Uptake, Wipe Tests, Schilling Test, Blood Volume, RBC Survival and Bio-Assay, and excellent technical support from our staff, you'll be up and running in no time.



**New Developments:** Capintec introduces the NumaCap™ software package; a DICOM/RIS software upgrade for new and existing Captus® 3000 Thyroid Uptake Systems. For a demonstration, visit our RSNA booth, or call us at (201) 825-9500 and ask a Capintec sales representative for more information!

### **Captus® 3000**

## Employee Updates

**Eric Offner** joined Capintec in November 1999 as an assembler in the CRC production area, becoming a technician shortly thereafter. In June 2000, Eric was promoted to Product Coordinator for the CRC dose calibrators. After three years, he was transferred to the TUS department as Product Coordinator. As Product Coordinator, Eric ensured that new orders were taken and repairs were made and returned to the customer in a timely fashion. These positions afforded Eric a well-rounded education of the equipment along with the opportunity to interact with customers by phone, email or at trade shows.

In August 2006, Eric joined the Product Support staff as Technical Sales & Applications Specialist. This new position allows Eric to contribute his extensive knowledge of the Captus and CRC line of products to our support team. In this new position, Eric will work directly with customers providing telephone applications as well as in the field with on-site installations and training. Eric will also provide support for product shows at regional and national meetings.


In September 2006, **Kathy Thomas**, MHA, CNMT, PET, joined the Capintec team as Regional Manager, Applications and Technical Support. Kathy comes to Capintec with 30+ years in nuclear medicine. Prior to joining Capintec, Kathy was a Senior Nuclear Medicine Technologist at City of Hope National Medical Center. She will continue to reside on the West coast to provide support to our customers in that region.

Kathy was awarded a Masters of Health Administration from LaVerne University and a BS from California Polytechnic State University. Her professional affiliations include the Society of Nuclear Medicine, American Healthcare Radiology Administrator, Nuclear Medicine Technology Certification Board, American Registry of Radiologic Technologists and the Technologist Advisory Board. Kathy has served as President of the Society of Nuclear Medicine Technologist Section.

Kathy's responsibilities will include providing management of technical and applications sales support to Capintec's West Coast customers as well as provide support for product shows at regional and national meetings ■

### Back by Popular Demand!

Capintec Periodic Table shirts:  
get one free at RSNA  
Booth #8929!



## Capintec Upcoming Show Schedule

Nov. 26 - Dec. 1  
**RSNA [Booth #8929]**  
Chicago, IL

### "In Step With Capintec" Staff

Layout & Editor.....Mark Nodzon  
Editor.....Margaret Corcoran

"In Step With Capintec" © 2006, Capintec, Inc. For your free subscription,  
write to: **In Step With Capintec, 6 Arrow Road, Ramsey, NJ 07446**

You can also view us on the web at [www.capintec.com](http://www.capintec.com)  
under our "News" section



# CAPINTEC, INC.

6 Arrow Road, Ramsey NJ 07446  
(201) 825-9500 or (800) 631-3826

## Lead-Lined Cabinetry from Capintec

Everything about Capintec's lead-lined cabinetry conveys quality. Attention to detail is what you expect from Capintec...and that's what you get. Manufactured using only the highest quality materials for a lifetime of reliable service. The use of heavy duty door hinges and drawer/shelf glides ensure effortless, trouble-free use. Its sleek, modern appearance with smooth, rounded corners and triple-coated paint finish will enhance any nuclear medicine department or PET facility. Capintec's lead-lined cabinetry has been

specifically designed for use in all areas where radionuclides and radioactive sources are stored, handled, or used. You can choose from a selection of standard cabinet configurations, or have lead-lined cabinets manufactured to your exact specifications. Members of our engineering and manufacturing teams are prepared to assist you in customizing lead-lined cabinetry to meet your requirements!

Visit Capintec booth #8929 at RSNA 2006  
for a free updated cabinetry brochure!

### Features:

- Lifetime guarantee on frame and assembly
- All units have solid steel frames, *an industry exclusive!*
- Specify your choice of lead shielding for all nuclear medicine and PET requirements (1/4", 1/2", 1"); select models available in 2" thickness for PET
- Easy glide...easy open & close draw, shelving and door design
- Drawers have a weight capacity of 100 lbs.
- Individual key locks for all drawers/doors
- Stainless steel countertops with spill-proof lip and backsplash
- Seamless single stainless steel countertop
- Rounded corners
- Free design service: custom cabinetry designed to your specifications

